



A.V Fistula Needles

Brief description of the product

The Fistula Needles are assembled by protective cap, needle tube, double-wing plate, lock fitting, tubing, inner conical interface, lock cover. All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free.

The Fistula Needles is intended to be used with blood composition collecting machines (for example centrifugalization style and rotating membrane style etc.) or blood dialysis machine for venous or arterial blood collecting work, then administer blood composition return to human body.

Characteristic

- ◆ Fixed type, rotating type, safety type
- ◆ Latex free
- ◆ Ultra-thin wall needles allow maximum blood flow
- ◆ Ultra-sharp for easier penetration
- ◆ Siliconized for smoother, more comfortable penetration
- ◆ Easy removal of luer cap
- ◆ Double entry points on needle for quick collection
- ◆ A.V. Fistula needles with safety device, can shield the cannula after use, preventing needlepoint from stabbing nurse, patient or other people, which can effectively reduce infection.

Detailed specifications

Type	SIZE		COLOR COLOR & CODE	LENGTH			ID RW	TW	REF
	OD(mm)	GAUGE		3/4"	1"	1 1/4"			
Fixed Type	1.40	17G	orange						01-901
	1.60	16G	green						01-902
	1.80	15G	beige						01-903
Rotating Type	1.40	17G	orange						01-904
	1.60	16G	green						01-905
	1.80	15G	beige						01-906
Safety Type	1.40	17G	orange						01-907
	1.60	16G	green						01-908
	1.80	15G	beige						01-909

Remark: normal size in above table, we can also make special size as clients' request.

Product Conformance



In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP,PC,PVC, SUS304 Stainless Steel Cannula,Silicone oil



Fixed Type

Rotating Type

Safety Type

