



## Irrigation Needles

### Brief description of the product

The Irrigation Needles are assembled by hub, needle tube, protect cap. All materials meet medical requirement; sterilized by ETO, pyrogen-free. The Irrigation Needles are intended for irrigation purposes in dentistry and ophthalmology, or used for the withdrawal of drugs/solutions out of vials and - if desired - subsequent injection into I.V. systems.



### Characteristic

- ◆ Multi angle bended, easy for cleanout, can be used for different occasion.
- ◆ Different type, client can choose products according to irrigation part.

### Detailed specifications

SIZE		COLOR	REF
OD (mm)	GAUGE	COLOR & CODE	
0.40	27G	Grey	01-1301
0.45	26G	brown	01-1302
0.50	25G	orange	01-1303
0.55	24G	purple	01-1304
0.60	23G	blue	01-1305
0.70	22G		01-1306
0.80	21G	green	01-1307
0.90	20G	yellow	01-1308
1.10	19G	cream	01-1309
1.20	18G	pink	01-1310

Remark: normal size in above table, we can also make special size as clients' request.



### Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: I)

### Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

### Main material

PP, SUS304 Stainless Steel Cannula, Silicone oil

