



Irrigation Syringes

Brief description of the product

The Irrigation Syringes are assembled by barrel, piston and plunger. All parts and material for this product meet medical requirements, after sterilized by ETO, pyrogen free.

The Irrigation Syringes are mainly used for hemorrhoids patient who need clean the anus before and after the operation, to speed up the resume of hurt place, fill in the catheter.

Characteristic

- ◆ Bulb type, ring type, flat type.
- ◆ Fully imported raw materials; barrel is transparent, easy to observe, scale printing ink adhesion is strong, not fall off.
- ◆ Edge generous, grip comfort, not easy to cause deformation.
- ◆ General cooperates: needle tube joint can match and gastric tube joint.

Detailed specifications

TYPE	REF
Bulb	02-701
Ring	02-702
Flat	02-703

Product Conformance

Conform to ISO 7886-1 and YZB/GUO 5665-66-2012

In compliance with European Medical Device Directive 93/42/EEC(CE Class: I)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, Isoprene rubber, Silicone oil, TPE

