



Safety Syringes

Brief description of the product

The Safety Syringes are assembled by protective cap, hypodermic needles, Shield shell, piston, barrel, plunger. All parts and material for this product meet medical requirements, after sterilized by ETO, pyrogen free.

The Safety Syringes used together with hypodermic needles are intended to inject the drug for the patient or take out the blood for the patient.

Characteristic

- ◆ Different size,meet different requirements.
- ◆ Central luer lock
- ◆ Special design to prevent plunger slipping out.
- ◆ Latex/latex free piston.
- ◆ W/needle,w/o needle.
- ◆ Custom needle size per specification.

Detailed specifications

TYPE	REF
1ml	02-801
3ml	02-802
5ml	02-803
10ml	02-804
20ml	02-805
50ml	02-806

Product Conformance

Conform to ISO 7886-1 and ISO 7886-4

In compliance with European Medical Device Directive 93/42/EEC(CE Class: W/needle: IIa;W/O needle: I)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, SUS304 Stainless Steel Cannula,Isoprene rubber,Silicone oil

