

# Chiba Needle

## Structure and composition:

The chiba needle consist of needle hub, needle tube Protective cap. All materials meet medical require ment, sterilized by ETO, pyrogen-free.

## Material:

PC,PP, SUS304 Stainless Steel Cannula,Silicone oil

The internal echogenic marker on the tip enables a correct placement of the needle and its constant visualization under ultrasound guidance.

Centimeter markings on the cannula surface enhance an easy determination of the insertion depth with the maximum safety for the patient.

Sharper tip design for easy puncture operation



The needle is coloured according to the international colour codes in order to help the user to identify the device gauge.

**Remark:** The type can make as clients' request.

## Intended use:

The Chiba needle iis used to inject drugs, necessary tools for lower guide wire and take liquid cell tissue fluid.

## Detailed specifications

Specifications	Length
18G	90mm、150mm、200mm
20G	90mm、150mm、200mm
22G	90mm、150mm、200mm

