


# Seldinger Needle

## Structure and composition:

The seldinger needle consist of needle hub, needle tube Protective cap. All materials meet medical require ment, sterilized by ETO, pyrogen-free.

## Material:

PC,PP, SUS304 Stainless Steel Cannula,Silicone oil




## Intended use:

The seldinger needle is used to penetrate into artery vas at the beginning of the intervention operation, and then introduce the guidewire into the vas.

## Detailed specifications

Specifications	Length
14G ~25G	90mm、150mm、200mm



**Remark:** The type can make as clients' request.