



Precise | Safety



Bone Marrow Biopsy Needles

SHANGHAI KINDLY(KDL) GROUP
NANCHANG KINDLY MEDICAL TECHNOLOGY CO.,LTD

F: +86-791 86861216 W: www.kdInc.com
T: +86-791 86861216

F/7, Tower A, No.77 North Aixihu Road,Nanchang High-Tech Industry Development Zone,Nanchang, Jiangxi 330096.China

The handle includes a sharpened luer-lock which gives a safer Attachment to the syringe. Aspiration can be performed at any Step of the procedure


Digital scale can provide accurate puncture depth

Multiple edge

The edge is strong and sharp, that allows a quick and less invasive penetration into the cortical

Special trocar needle and sharp tubing ensure minimal invasion, 5 side holes at front end of tubing, multi-way collect bone marrow, more efficient and less pain.

Pick up needle is provided with a rotation coupling security system which limits the risk of movement of the stylet during the insertion.



Structure and composition:

The Bone marrow biopsy needle consist of needle hub, needle tube, lining core, handle and Protective cap.

Material:

ABS,PP PC,SUS304 Stainless Steel Cannula,Silicone oil

Intended use:

The Bone marrow biopsy needle is used for bone marrow aspiration.



Specification

SIZE		LENGTH			
OD(mm)	GAUGE	70 2 ¾"	90 3 ½"	110 3 ¼"	130 5"
1.8	15G	✓		✓	
2.5	13G	✓		✓	✓
3.0	11G		✓		

Remark: normal size in above table, we can also make special size as clients' request.