




Huber Needles

Brief description of the product

The Huber Needles are assembled by protective cap, needle hub, needle tube, tubing, injection site, robert clamp and other parts. All materials meet medical requirement, sterilized by ETO, pyrogen-free. It's intended to be used to feed medicine to the device implanted in patient hypoderm.

Characteristic

- ◆ Latex free
- ◆ Non coring
- ◆ 90 degrees SUS304 stainless steel needle
- ◆ A rotating wing provides less stress on both the tissue and the port
- ◆ Occlusion clamp



Detailed specifications

SIZE		COLOR	LENGTH								ID		Bevel angle		REF
OD (mm)	GAUGE	COLOR & CODE	1/3"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2"	RW	TW	Long bevel	Short bevel	
0.70	22G	black													01-1001
0.80	21G	green													01-1002
0.90	20G	yellow													01-1003
1.10	19G	cream													01-1004

Remark: normal size in above table, we can also make special size as clients' request.

Product Conformance


In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PC,PP,PVC, SUS304 Stainless Steel Cannula,Silicone oil



1 Scalp Vein Sets Type

2 Hypodermic Needles Type

