



Anaesthesia Needles

Extradural Anaesthesia Needles(AN-E)

- ◆ Bend and round-smooth needlepoint can greatly reduce the risks of breaking the hard spinal film and assure the cannula entering successfully.
- ◆ Graduation on cannula, convenient for control.
- ◆ Needle hub with wings, convenient to be hold and needle sticking into easily.



Spinal Needles(AN-SI)

- ◆ Transparent hub, easy to observe the cerebrospinal fluid discharge.
- ◆ Hub with lock can keep the surface of needle bevel joint well.



Non-hurt Spinal Needles(AN-SII)

- ◆ Special design for needlepoint will not hurt the hard spinal film, close puncture hole automatically and reduce cerebrospinal fluid discharge.
- ◆ Suitable shape and location of side hole ensure enough flux, meanwhile reduce the happening of faulty retardance.




Nerve Anaesthesia Needles(AN-N)

- ◆ Used for part anaesthesia, ease pain or injecting local anesthetic when in emergency.
- ◆ For single use only, avoid cross infection.



Nerve Anaesthesia Needles(AN-X)


- ◆ Bend and round-smooth needlepoint can greatly reduce the risks of breaking the hard spinal film and assure the cannula entering successfully.
- ◆ Graduation on cannula, convenient for control.
- ◆ Needle hub with wings, convenient to be hold and needle sticking into easily.
- ◆ Small OD, greatly reduce the hurt on pediatric patients.
- ◆ Soft tube connecting with needle hub, Easily use and safety.



Detailed specifications

Type	SIZE		COLOR COLOR&CODE	LENGTH						ID		REF
	OD(mm)	GAUGE		2 1/4	2 1/2	3 1/8	3 1/2	4 1/2	4	RW	TW	
AN-E Extradural Anaesthesia Needles	0.80	21G	green									01-701
	0.90	20G	yellow									01-702
	1.00	19G	cream									01-703
	1.20	18G	pink									01-704
	1.60	16G	white		white	white	white				white	01-705
	1.80	15G	Blue-grey									01-706
AN-S I Spinal Needles	0.40	27G	Grey									01-707
	0.45	26G	brown									01-708
	0.50	25G	orange									01-709
	0.55	24G	purple									01-710
	0.60	23G	blue									01-711
	0.70	22G	black									01-712
	0.80	21G	green									01-713
	0.90	20G	yellow									01-714
	1.00	19G	cream									01-715
	1.20	18G	pink									01-716
AN-S II No-hurt Spinal Needles	0.50	25G	orange									01-717
	0.55	24G	purple									01-718
	0.60	23G	blue									01-719
	0.70	22G	black									01-720
	0.80	21G	green									01-721
	0.90	20G	yellow									01-722
	1.10	19G	cream									01-723
	1.20	18G	pink									01-724
AN-N Nerve Anaesthesia Needles	0.50	25G	orange									01-725
	0.55	24G	purple									01-726
	0.60	23G	blue									01-727
	0.70	22G	black									01-728
	0.80	21G	green									01-729
	0.90	20G	yellow									01-730
AN-X Pediatric Anesthesia Needles	0.50	25G	orange	3/4								01-731
	0.55	24G	purple	3/4								01-732
												01-733

Remark: normal size in above table, we can also make special size as clients' request.



1 (AN-E)

2 (AN-SI)

3 (AN-SII)

4 (AN-N)

5 (AN-X)

Product Conformance

Conform to BS6196

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa,III)

Main material

PP,PC,ABS, SUS304 Stainless Steel Cannula,Silicone oil

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.