

# Pre-Filled Flush Syringes



## Improve Clinical Workflow and Productivity by:

- Reducing the number of steps and time required for clinicians to prepare flush syringes.
- Reducing acquisition, storage and distribution costs for the multiple components currently required to flush IV sites.
- Reducing sharps collector waste by up to 50%, compared to a conventional 12ml syringe.

## Cost-Effective

- Eliminate flush syringe filling labour
- Reduce costs for drug wastage
- Simplification of purchasing and inventory control



Unfilled brine



TYPE	REF
3ml	02-1001
5ml	02-1002
10ml	02-1003

The Pre-Filled Flush Syringes are assembled by barrel, plunger, piston, syringe cap and Saline. All parts and material for this product meet medical requirements, after sterilized by moist heat. It's used to seal and wash the catheter end during the internal of different drug treatment.

Single-use pre-filled syringes help reduce medication errors and contamination risks associated with multi-dose vials. It also improves clinician efficiency and workflow by eliminating the steps and time involved in the manual preparation of saline syringes.

## Reduce Medical Errors and Complications by:

- Eliminating potential cross-contamination associated with multi-dose containers.
- Decreasing breaks in aseptic technique during preparation.
- Using preservative-free saline, suitable for use with all patients.



## Convenient

- Less preparation, fewer materials, easy storage and disposal
- Remove need for vials, ampoules and needles with a single product
- Simplify catheter maintenance procedures

## Characteristic

- Three sizes comply different flush application
- BIIR piston, latex free
- Special design to prevent plunger slipping out.

**Luer lock**  
Firm and not easy to fall off



SHANGHAI KINDLY(KDL) GROUP  
NANCHANG KINDLY MEDICAL TECHNOLOGY CO.LTD

F7, Tower A, No.77 North Aixihu Road,Nanchang  
High-Tech Industry Development Zone,Nanchang,  
Jiangxi 330096,China  
F: +86-791 86861216  
T: +86-791 86861216  
W: www.kdlinc.com

