



Safety Scalp Vein Sets

Brief description of the product

The Safety Scalp Vein Sets are assembled by protective cap, needle tube, needle handle, tubing, needle hub, safety devices. All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free.

The Safety Scalp Vein Sets intended to be used with Syringes or Infusion/Transfusion sets for medicine infusion and blood collection.


Characteristic

- ◆ Safety type
- ◆ Soft wings
- ◆ Thin wall cannula, large ID, and smooth flow
- ◆ Safety Scalp Vein Sets, additional features
- ◆ Color coded wings by size for clear recognition
- ◆ Protective sheath on the tube
- ◆ Easy activation of safety sheath, prevent accidental needle sticks
- ◆ Needle hub: 1. PVC hub(slip with cap)
2. ABS hub(luer lock with cap)
- ◆ Packing: in bulk, PE bag, in blister.

Detailed specifications

SIZE	COLOR	LENGTH				ID		Bevel angle		REF		
		OD (mm)	GAUGE	COLOR & CODE	1/2	5/8	3/4	1	RW		TW	Long bevel
0.40	27G	Grey										01-301
0.45	26G	brown										01-302
0.50	25G	orange										01-303
0.55	24G	purple										01-304
0.60	23G	blue										01-305
0.70	22G	black										01-306
0.80	21G	green										01-307
0.90	20G	yellow										01-308
1.10	19G	cream										01-309
1.20	18G	pink										01-310

Remark: normal size in above table, we can also make special size as clients' request.



Comparison between Regular Wall, Thin Wall and Extra Thin Wall



Long bevel



Short bevel



Product Conformance

Conform to ISO 7864 and ISO 23908

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, SUS304 Stainless Steel Cannula, Silicone oil

